

EBP07503 - IMMUNIZING PEPTIDE

SKU: EBP07503

SPECIFICATIONS

Formulation	100µg of Lyophilised peptide, formulation De-ionized (DI) water.
Unit Size	100 µg
Storage	
Instructions	Shipped at ambient temperature, store at -20°C.
Accession ID	NP_000168.1; NP_937895.1; NP_001121138.1; NP_001244958.1; NP_001244959.1; NP_001121135.2; NP_001340005.1; NP_001340007.1
Immunogen	Peptide with sequence C-PRLKDKKMDAHP, from the internal region of the protein sequence according to NP_000168.1; NP_937895.1; NP_001121138.1; NP_001244958.1; NP_001244959.1; NP_001121135.2; NP_001340005.1; NP_001340007.1.
Product Comments	The peptide was used in the production of EB07503 - Goat Anti-Gelsolin / GSN Antibody. It is currently untested in blocking applications.
Peptide Sequence	C-PRLKDKKMDAHP
Shipping Instructions	Ambient
Product Grade	https://prod-vector-labs-pimcore-assets.s3.us-east-1.amazonaws.com/assets/products/image/elite_medium.png

DOCUMENTS

- [Data Sheet](#)