

Email: <a href="mailto:customerservice@vectorlabs.com">customerservice@vectorlabs.com</a>

Telephone: (650) 697-3600

## **EBP07503 - IMMUNIZING PEPTIDE**

**SKU:** EBP07503

## **SPECIFICATIONS**

 $\textbf{Formulation} \ 100 \mu g \ of \ Lyophilised \ peptide, \ formulation \ De-ionized \ (DI) \ water.$ 

**Unit Size** 100 μα

Shipped at ambient temperature, store at -20°C.

 $NP\_000168.1; NP\_937895.1; NP\_001121138.1; NP\_001244958.1; NP\_001244959.1; NP\_001121135.2; NP\_001244959.1; NP\_001121135.2; NP\_00112115.2; NP\_00112115.2; NP\_00112115.2; NP\_00112115.2; NP\_00112115.2; NP\_00112115.2; NP\_00112115.2; NP\_001121$ Accession

NP\_001340005.1; NP\_001340007.1 ID

Peptide with sequence C-PRLKDKKMDAHP, from the internal region of the protein sequence according to

Immunogen NP 000168.1; NP 937895.1; NP 001121138.1; NP 001244958.1; NP 001244959.1; NP 001121135.2;

NP 001340005.1; NP 001340007.1.

The peptide was used in the production of EB07503 - Goat Anti-Gelsolin / GSN Antibody. It is currently untested **Product** 

**Comments** in blocking applications.

**Peptide** 

C-PRLKDKKMDAHP Sequence

Shipping

Instructions Ambient

**Product** 

https://prod-vector-labs-pimcore-assets.s3.us-east-1.amazonaws.com/assets/products/image/elite\_medium.png Grade