

Email: customerservice@vectorlabs.com

Telephone: (650) 697-3600

EBP07515 - IMMUNIZING PEPTIDE

SKU: EBP07515

SPECIFICATIONS

Formulation 100µg of Lyophilised peptide, formulation De-ionized (DI) water.

Unit Size 100 μg

Storage

Instructions Shipped at ambient temperature, store at -20°C.

Accession

NP_034311.1 ID

Peptide with sequence C-KELRKLEEHNEDD, from the internal region (near the C Terminus) of the protein sequence

according to NP_034311.1.

The peptide was used in the production of EB07515 - Goat Anti-Fibrillin 2 (mouse) Antibody. It is currently Product

Comments untested in blocking applications.

Peptide

C-KELRKLEEHNEDD Sequence

Shipping

Ambient Instructions

Product

https://prod-vector-labs-pimcore-assets.s3.us-east-1.amazonaws.com/assets/products/image/aspiring_medium.png Grade