

EBP07520 - IMMUNIZING PEPTIDE

SKU: EBP07520



SPECIFICATIONS

Formulation	100µg of Lyophilised peptide, formulation De-ionized (DI) water.
Unit Size	100 µg
Storage Instructions	Shipped at ambient temperature, store at -20°C.
Accession ID	NP_055306.1; NP_683696.2; NP_683697.2; NP_683698.2; NP_001166237.1
Immunogen	Peptide with sequence C-DEVEYQKRRSQKIT, from the internal region of the protein sequence according to NP_055306.1; NP_683696.2; NP_683697.2; NP_683698.2; NP_001166237.1.
Product Comments	The peptide was used in the production of EBP07520 - Goat Anti-FOXP2 (internal) Antibody. It is currently untested in blocking applications.
Peptide Sequence	C-DEVEYQKRRSQKIT
Shipping Instructions	Ambient
Product Grade	https://prod-vector-labs-pimcore-assets.s3.us-east-1.amazonaws.com/assets/products/image/elite_medium.png

DOCUMENTS

- [Data Sheet](#)