

Email: <a href="mailto:customerservice@vectorlabs.com">customerservice@vectorlabs.com</a>

Telephone: (650) 697-3600

## **EBP07540 - IMMUNIZING PEPTIDE**

**SKU:** EBP07540

## **SPECIFICATIONS**

 $\textbf{Formulation} \ 100 \mu g \ of \ Lyophilised \ peptide, \ formulation \ De-ionized \ (DI) \ water.$ 

**Unit Size** 100 μα

Storage

Shipped at ambient temperature, store at -20°C. Instructions

**Accession** 

ID

Peptide with sequence KSCREKTELTPDQQ, from the internal region of the protein sequence according to

Immunogen NP 005114.1; NP 001193906.1; NP 001193922.1; NP 001193921.1; NP 001193907.1.

NP 005114.1; NP 001193906.1; NP 001193922.1; NP 001193921.1; NP 001193907.1

The peptide was used in the production of EB07540 - Goat Anti-Farnesoid X receptor / FXR Antibody. It is **Product** 

**Comments** currently untested in blocking applications.

**Peptide** 

**KSCREKTELTPDQQ** Sequence

Shipping

Instructions Ambient

**Product** 

Grade

https://prod-vector-labs-pimcore-assets.s3.us-east-1.amazonaws.com/assets/products/image/elite\_medium.png