



# EBP07587-B - IMMUNIZING PEPTIDE

**SKU:** EBP07587-B



---

## SPECIFICATIONS

|                              |   |
|------------------------------|---|
| <b>Formulation</b>           | 100µg of Lyophilised peptide, formulation De-ionized (DI) water.  |
| <b>Unit Size</b>             | 100 µg  |
| <b>Storage Instructions</b>  | Shipped at ambient temperature, store at -20°C.   |
| <b>Accession ID</b>          | NP_003191.1   |
| <b>Immunogen</b>             | Peptide with sequence C-KAPRARTSPDEDED., from the internal region of the protein sequence according to NP_003191.1.   |
| <b>Product Comments</b>      | The peptide was used in the production of EBP07587-B - Goat Anti-TCF3 / ITF1, Biotinylated Antibody. It is currently untested in blocking applications.   |
| <b>Peptide Sequence</b>      | C-KAPRARTSPDEDED.   |
| <b>Shipping Instructions</b> | Ambient   |
| <b>Product Grade</b>         | <a href="https://prod-vector-labs-pimcore-assets.s3.us-east-1.amazonaws.com/assets/products/image/elite_medium.png">https://prod-vector-labs-pimcore-assets.s3.us-east-1.amazonaws.com/assets/products/image/elite_medium.png</a> |

## DOCUMENTS

- [Data Sheet](#)