

Email: customerservice@vectorlabs.com

Telephone: (650) 697-3600

EBP07595 - IMMUNIZING PEPTIDE

SKU: EBP07595

SPECIFICATIONS

Formulation 100µg of Lyophilised peptide, formulation De-ionized (DI) water.

Unit Size 100 μg

Storage

Instructions Shipped at ambient temperature, store at -20°C.

Accession

NP_061834.1; NP_001001895.1 ID

 $\label{eq:local_$

The peptide was used in the production of EB07595 - Goat Anti-STS2 / TULA (N Terminus) Antibody. It is currently Product

Comments untested in blocking applications.

Peptide

QLYAKVSNKLKSRS-C Sequence

Shipping

Ambient Instructions

Product

https://prod-vector-labs-pimcore-assets.s3.us-east-1.amazonaws.com/assets/products/image/aspiring_medium.png Grade