

EBP07600 - IMMUNIZING PEPTIDE

SKU: EBP07600



SPECIFICATIONS

Formulation	100µg of Lyophilised peptide, formulation De-ionized (DI) water.
Unit Size	100 µg
Storage	Shipped at ambient temperature, store at -20°C.
Instructions	
Accession ID	NP_055099.1; NP_001014841.1
Immunogen	Peptide with sequence C-HYLLQVGSEKQKEP, from the internal region of the protein sequence according to NP_055099.1; NP_001014841.1.
Product Comments	The peptide was used in the production of EB07600 - Goat Anti-Neurochondrin Antibody. It is currently untested in blocking applications.
Peptide Sequence	C-HYLLQVGSEKQKEP
Shipping Instructions	Ambient
Product Grade	https://prod-vector-labs-pimcore-assets.s3.us-east-1.amazonaws.com/assets/products/image/elite_medium.png

DOCUMENTS

- [Data Sheet](#)