

EBP07646 - IMMUNIZING PEPTIDE

SKU: EBP07646



SPECIFICATIONS

Formulation 100µg of Lyophilised peptide, formulation De-ionized (DI) water.

Unit Size 100 µg

Storage Instructions Shipped at ambient temperature, store at -20°C.

Accession ID NP_085048.1

Immunogen Peptide with sequence C-NYAKVRKMPET, from the internal region of the protein sequence according to NP_085048.1.

Product Comments The peptide was used in the production of EBP07646 - Goat Anti-NDFIP1 Antibody. It is currently untested in blocking applications.

Peptide Sequence C-NYAKVRKMPET

Shipping Instructions Ambient

Product Grade https://prod-vector-labs-pimcore-assets.s3.us-east-1.amazonaws.com/assets/products/image/aspiring_medium.png