

EBP07649 - IMMUNIZING PEPTIDE

SKU: EBP07649



SPECIFICATIONS

Formulation	100µg of Lyophilised peptide, formulation De-ionized (DI) water.
Unit Size	100 µg
Storage	Shipped at ambient temperature, store at -20°C.
Instructions	
Accession ID	NP_849180.1; NP_000448.3; NP_787110.2; NP_001025174.1; NP_001245284.1; NP_001274111.1; NP_001274112.1
Immunogen	Peptide with sequence C-HPHLMQEHEMGTN, from the internal region of the protein sequence according to NP_849180.1; NP_000448.3; NP_787110.2; NP_001025174.1; NP_001245284.1; NP_001274111.1; NP_001274112.1.
Product Comments	The peptide was used in the production of EBP07649 - Goat Anti-HNF4A Antibody. It is currently untested in blocking applications.
Peptide Sequence	C-HPHLMQEHEMGTN
Shipping Instructions	Ambient
Product Grade	https://prod-vector-labs-pimcore-assets.s3.us-east-1.amazonaws.com/assets/products/image/elite_medium.png

DOCUMENTS

- [Data Sheet](#)