

EBP07687 - IMMUNIZING PEPTIDE

SKU: EBP07687



SPECIFICATIONS

| | |
|------------------------------|--|
| Formulation | 100µg of Lyophilised peptide, formulation De-ionized (DI) water. |
| Unit Size | 100 µg |
| Storage | |
| Instructions | Shipped at ambient temperature, store at -20°C. |
| Accession ID | NP_006425.2; NP_056200.1; NP_001030126.1; NP_001030127.1; NP_001030128.1; NP_079267.1; NP_001030129.1; NP_001277225.2; NP_001364131.1; NP_001371378.1 |
| Immunogen | Peptide with sequence C-SERNDWEPPDKKVD, from the internal region of the protein sequence according to NP_006425.2; NP_056200.1; NP_001030126.1; NP_001030127.1; NP_001030128.1; NP_079267.1; NP_001030129.1; NP_001277225.2; NP_001364131.1; NP_001371378.1. |
| Product Comments | The peptide was used in the production of EB07687 - Goat Anti-SORBS1/ponsin (internal region) Antibody. It is currently untested in blocking applications. |
| Peptide Sequence | C-SERNDWEPPDKKVD |
| Shipping Instructions | Ambient |
| Product Grade | https://prod-vector-labs-pimcore-assets.s3.us-east-1.amazonaws.com/assets/products/image/elite_medium.png |

DOCUMENTS

- [Data Sheet](#)