

EBP07687 - IMMUNIZING PEPTIDE

SKU: EBP07687



SPECIFICATIONS

Formulation	100µg of Lyophilised peptide, formulation De-ionized (DI) water.
Unit Size	100 µg
Storage	
Instructions	Shipped at ambient temperature, store at -20°C.
Accession ID	NP_006425.2; NP_056200.1; NP_001030126.1; NP_001030127.1; NP_001030128.1; NP_079267.1; NP_001030129.1; NP_001277225.2; NP_001364131.1; NP_001371378.1
Immunogen	Peptide with sequence C-SERNDWEPPDKKVD, from the internal region of the protein sequence according to NP_006425.2; NP_056200.1; NP_001030126.1; NP_001030127.1; NP_001030128.1; NP_079267.1; NP_001030129.1; NP_001277225.2; NP_001364131.1; NP_001371378.1.
Product Comments	The peptide was used in the production of EBP07687 - Goat Anti-SORBS1/ponsin (internal region) Antibody. It is currently untested in blocking applications.
Peptide Sequence	C-SERNDWEPPDKKVD
Shipping Instructions	Ambient
Product Grade	https://prod-vector-labs-pimcore-assets.s3.us-east-1.amazonaws.com/assets/products/image/elite_medium.png

DOCUMENTS

- [Data Sheet](#)