

Email: customerservice@vectorlabs.com

Telephone: (650) 697-3600

## EBP07687 - IMMUNIZING PEPTIDE

**SKU:** EBP07687

## **SPECIFICATIONS**

 $\textbf{Formulation} \ 100 \mu g \ of \ Lyophilised \ peptide, \ formulation \ De-ionized \ (DI) \ water.$ 

**Unit Size** 100 μα

Shipped at ambient temperature, store at -20°C.

NP\_006425.2; NP\_056200.1; NP\_001030126.1; NP\_001030127.1; NP\_001030128.1; NP\_079267.1; Accession

NP\_001030129.1; NP\_001277225.2; NP\_001364131.1; NP\_001371378.1 ID

Peptide with sequence C-SERNDWEPPDKKVD, from the internal region of the protein sequence according to

Immunogen NP 006425.2; NP 056200.1; NP 001030126.1; NP 001030127.1; NP 001030128.1; NP 079267.1;

NP\_001030129.1; NP\_001277225.2; NP\_001364131.1; NP\_001371378.1.

The peptide was used in the production of EB07687 - Goat Anti-SORBS1/ponsin (internal region) Antibody. It is **Product** 

**Comments** currently untested in blocking applications.

**Peptide** 

C-SERNDWEPPDKKVD Sequence

Shipping

Instructions Ambient

**Product** 

https://prod-vector-labs-pimcore-assets.s3.us-east-1.amazonaws.com/assets/products/image/elite\_medium.png Grade