

EBP07691 - IMMUNIZING PEPTIDE

SKU: EBP07691



SPECIFICATIONS

Formulation	100µg of Lyophilised peptide, formulation De-ionized (DI) water.
Unit Size	100 µg
Storage Instructions	Shipped at ambient temperature, store at -20°C.
Accession ID	NP_005154.2; NP_001014431.1; NP_001014432.1
Immunogen	Peptide with sequence C-QDVDQREAPLN, from the internal region of the protein sequence according to NP_005154.2; NP_001014431.1; NP_001014432.1.
Product Comments	The peptide was used in the production of EB07691 - Goat Anti-AKT1 Antibody. It is currently untested in blocking applications.
Peptide Sequence	C-QDVDQREAPLN
Shipping Instructions	Ambient
Product Grade	https://prod-vector-labs-pimcore-assets.s3.us-east-1.amazonaws.com/assets/products/image/elite_plus_medium.png

DOCUMENTS

- [Data Sheet](#)