

Email: customerservice@vectorlabs.com

Telephone: (650) 697-3600

EBP07697 - IMMUNIZING PEPTIDE

SKU: EBP07697

SPECIFICATIONS

Formulation 100µg of Lyophilised peptide, formulation De-ionized (DI) water.

Unit Size 100 μg

Storage Shipped at ambient temperature, store at -20°C.

Accession NP_000041.2; NP_446464.1 ID

Immunogen Peptide with sequence C-ENPKNQAPPGLYTKTQD, from the internal region of the protein sequence according to NP_000041.2; NP_446464.1.

The peptide was used in the production of EB07697 - Goat Anti-Argininosuccinate synthetase 1 Antibody. It is **Product**

Comments currently untested in blocking applications.

Peptide

C-ENPKNQAPPGLYTKTQD Sequence

Shipping

Ambient Instructions

Product

 $https://prod-vector-labs-pimcore-assets.s3.us-east-1.amazonaws.com/assets/products/image/elite_plus_medium.png$ Grade