

Email: customerservice@vectorlabs.com

Telephone: (650) 697-3600

EBP07706 - IMMUNIZING PEPTIDE

SKU: EBP07706

SPECIFICATIONS

Formulation 100µg of Lyophilised peptide, formulation De-ionized (DI) water.

Unit Size 100 μg

Storage Shipped at ambient temperature, store at -20°C.

Accession

NP_000036.2; NP_001231367.1; NP_001355949.1 ID

Peptide with sequence C-REGNHKPIDYLNPPK, from the C Terminus of the protein sequence according to

Immunogen NP_000036.2; NP_001231367.1; NP_001355949.1.

The peptide was used in the production of EB07706 - Goat Anti-Arginase, type 1 / ARG1 Antibody. It is currently **Product**

Comments untested in blocking applications.

Peptide

C-REGNHKPIDYLNPPK Sequence

Shipping

Ambient Instructions

Product

 $https://prod-vector-labs-pimcore-assets.s 3. us-east-1. amazonaws. com/assets/products/image/elite_plus_medium.png$ Grade