

Email: <a href="mailto:customerservice@vectorlabs.com">customerservice@vectorlabs.com</a>

Telephone: (650) 697-3600

## **EBP07716 - IMMUNIZING PEPTIDE**

**SKU:** EBP07716

## **SPECIFICATIONS**

 $\textbf{Formulation} \ 100 \mu g \ of \ Lyophilised \ peptide, \ formulation \ De-ionized \ (DI) \ water.$ 

100 μg **Unit Size** 

Storage

Shipped at ambient temperature, store at -20°C. Instructions

**Accession** 

NP 003346.2 ID

Peptide with sequence C-DSVKQFYTKGSEH, from the internal region of the protein sequence according to **Immunogen** 

NP\_003346.2.

The peptide was used in the production of EB07716 - Goat Anti-Uncoupling protein 2 / UCP2 Antibody. It is **Product** 

**Comments** currently untested in blocking applications.

**Peptide** 

C-DSVKQFYTKGSEH Sequence

**Shipping** 

Instructions Ambient

**Product** 

https://prod-vector-labs-pimcore-assets.s3.us-east-1.amazonaws.com/assets/products/image/elite\_medium.png Grade