

Email: <a href="mailto:customerservice@vectorlabs.com">customerservice@vectorlabs.com</a>

Telephone: (650) 697-3600

## **EBP07721 - IMMUNIZING PEPTIDE**

**SKU:** EBP07721

## **SPECIFICATIONS**

Formulation 100µg of Lyophilised peptide, formulation De-ionized (DI) water.

**Unit Size** 100 μg

Storage

Instructions Shipped at ambient temperature, store at -20°C.

Accession

NP\_033363.1; NP\_001276397.1 ID

 $\label{eq:local_$ 

The peptide was used in the production of EB07721 - Goat Anti-Tcl1 (mouse) Antibody. It is currently untested in **Product** 

**Comments** blocking applications.

**Peptide** 

Sequence

**QLYPKDRYRSCD** 

Shipping

Ambient Instructions

**Product** 

Grade

https://prod-vector-labs-pimcore-assets.s3.us-east-1.amazonaws.com/assets/products/image/aspiring\_medium.png