

EBP07727 - IMMUNIZING PEPTIDE

SKU: EBP07727



SPECIFICATIONS

Formulation	100µg of Lyophilised peptide, formulation De-ionized (DI) water.
Unit Size	100 µg
Storage	
Instructions	Shipped at ambient temperature, store at -20°C.
Accession ID	NP_002657.3; NP_001138783.1 NP_001138783.1
Immunogen	Peptide with sequence C-PREKPKRRVSDS, from the internal region (near the C Terminus) of the protein sequence according to NP_002657.3; NP_001138783.1 NP_001138783.1.
Product	The peptide was used in the production of EBP07727 - Goat Anti-Perilipin 1 (internal) Antibody. It is currently
Comments	untested in blocking applications.
Peptide Sequence	C-PREKPKRRVSDS
Shipping Instructions	Ambient
Product Grade	https://prod-vector-labs-pimcore-assets.s3.us-east-1.amazonaws.com/assets/products/image/elite_medium.png

DOCUMENTS

- [Data Sheet](#)