

EBP07734 - IMMUNIZING PEPTIDE

SKU: EBP07734



SPECIFICATIONS

Formulation	100µg of Lyophilised peptide, formulation De-ionized (DI) water.
Unit Size	100 µg
Storage	Shipped at ambient temperature, store at -20°C.
Instructions	
Accession ID	NP_005929.2
Immunogen	Peptide with sequence C-EAASQDRMPQDLD, from the internal region of the protein sequence according to NP_005929.2.
Product Comments	The peptide was used in the production of EBP07734 - Goat Anti-FOXO4 / MLLT7 Antibody. It is currently untested in blocking applications.
Peptide Sequence	C-EAASQDRMPQDLD
Shipping Instructions	Ambient
Product Grade	https://prod-vector-labs-pimcore-assets.s3.us-east-1.amazonaws.com/assets/products/image/elite_medium.png

DOCUMENTS

- [Data Sheet](#)