

Email: customerservice@vectorlabs.com

Telephone: (650) 697-3600

EBP07737 - IMMUNIZING PEPTIDE

SKU: EBP07737

SPECIFICATIONS

 $\textbf{Formulation} \ 100 \mu g \ of \ Lyophilised \ peptide, \ formulation \ De-ionized \ (DI) \ water.$

Unit Size 100 μα

Storage

Shipped at ambient temperature, store at -20°C. Instructions

Accession

NP 067021.1; NP 872580.1 ID

Peptide with sequence C-RGKGDQKSAASQK, from the internal region of the protein sequence according to Immunogen

NP 067021.1; NP 872580.1.

The peptide was used in the production of EB07737 - Goat Anti-CTDSP1 Antibody. It is currently untested in **Product**

Comments blocking applications. **Peptide**

Sequence

C-RGKGDQKSAASQK

Shipping

Instructions Ambient

Product

https://prod-vector-labs-pimcore-assets.s3.us-east-1.amazonaws.com/assets/products/image/elite_medium.png Grade