

EBP07743 - IMMUNIZING PEPTIDE

SKU: EBP07743



SPECIFICATIONS

Formulation	100µg of Lyophilised peptide, formulation De-ionized (DI) water.
Unit Size	100 µg
Storage	Shipped at ambient temperature, store at -20°C.
Instructions	
Accession ID	NP_004466.2
Immunogen	Peptide with sequence C-QIYQDRDQFAK, from the internal region of the protein sequence according to NP_004466.2.
Product	The peptide was used in the production of EBP07743 - Goat Anti-Flotillin 2 / FLOT2 (internal) Antibody. It is currently untested in blocking applications.
Comments	
Peptide Sequence	C-QIYQDRDQFAK
Shipping Instructions	Ambient
Product Grade	https://prod-vector-labs-pimcore-assets.s3.us-east-1.amazonaws.com/assets/products/image/elite_medium.png