

Email: customerservice@vectorlabs.com

Telephone: (650) 697-3600

EBP07755 - IMMUNIZING PEPTIDE

SKU: EBP07755

SPECIFICATIONS

 $\textbf{Formulation} \ 100 \mu g \ of \ Lyophilised \ peptide, \ formulation \ De-ionized \ (DI) \ water.$

Unit Size 100 μα

Storage

Shipped at ambient temperature, store at -20°C. Instructions

Accession

NP 002582.2 ID

Peptide with sequence C-EKEVEDIEKYLEDQ, from the internal region (near the C Terminus) of the protein **Immunogen**

sequence according to NP 002582.2.

The peptide was used in the production of EB07755 - Goat Anti-PCK1 / PEPCKC (internal region) Antibody. It is **Product**

Comments currently untested in blocking applications.

Peptide

C-EKEVEDIEKYLEDQ Sequence

Shipping

Instructions Ambient

Product

https://prod-vector-labs-pimcore-assets.s3.us-east-1.amazonaws.com/assets/products/image/elite_medium.png Grade