

Email: <a href="mailto:customerservice@vectorlabs.com">customerservice@vectorlabs.com</a>

Telephone: (650) 697-3600

## **EBP07756 - IMMUNIZING PEPTIDE**

**SKU:** EBP07756

## **SPECIFICATIONS**

Formulation 100µg of Lyophilised peptide, formulation De-ionized (DI) water.

**Unit Size** 100 μg

Storage

Instructions Shipped at ambient temperature, store at -20°C.

Accession

NP\_149097.1 ID

Immunogen Peptide with sequence C-HLPQMGDQDGKRE, from the internal region of the protein sequence according to NP\_149097.1.

The peptide was used in the production of EB07756 - Goat Anti-GALP Antibody. It is currently untested in blocking **Product** 

Comments applications.

Peptide

C-HLPQMGDQDGKRE Sequence

Shipping

Ambient Instructions

**Product** 

https://prod-vector-labs-pimcore-assets.s3.us-east-1.amazonaws.com/assets/products/image/aspiring\_medium.png Grade