

Email: customerservice@vectorlabs.com

Telephone: (650) 697-3600

EBP07771 - IMMUNIZING PEPTIDE

SKU: EBP07771

SPECIFICATIONS

 $\textbf{Formulation} \ 100 \mu g \ of \ Lyophilised \ peptide, \ formulation \ De-ionized \ (DI) \ water.$

Unit Size 100 μα

Storage

Shipped at ambient temperature, store at -20°C. Instructions

Accession

NP 001049.1*; NP 056542.1* ID

Peptide with sequence C-EIPGDSSDRYHEQ, from the internal region of the protein sequence according to Immunogen

NP 001049.1*; NP 056542.1*.

The peptide was used in the production of EB07771 - Goat Anti-TACR1 Antibody. It is currently untested in **Product**

Comments blocking applications.

Peptide C-EIPGDSSDRYHEQ

Sequence

Shipping Instructions Ambient

Product Grade

https://prod-vector-labs-pimcore-assets.s3.us-east-1.amazonaws.com/assets/products/image/elite_medium.png