

EBP07782 - IMMUNIZING PEPTIDE

SKU: EBP07782



SPECIFICATIONS

Formulation	100µg of Lyophilised peptide, formulation De-ionized (DI) water.
Unit Size	100 µg
Storage Instructions	Shipped at ambient temperature, store at -20°C.
Accession ID	NP_780512.2; NP_958797.2; NP_958798.2
Immunogen	Peptide with sequence C-DQSKKGKGKPQN, from the internal region of the protein sequence according to NP_780512.2; NP_958797.2; NP_958798.2.
Product Comments	The peptide was used in the production of EBP07782 - Goat Anti-Sall4 (mouse) Antibody. It is currently untested in blocking applications.
Peptide Sequence	C-DQSKKGKGKPQN
Shipping Instructions	Ambient
Product Grade	https://prod-vector-labs-pimcore-assets.s3.us-east-1.amazonaws.com/assets/products/image/elite_medium.png

DOCUMENTS

- [Data Sheet](#)