

Email: customerservice@vectorlabs.com

Telephone: (650) 697-3600

EBP07782 - IMMUNIZING PEPTIDE

SKU: EBP07782

SPECIFICATIONS

 $\textbf{Formulation} \ 100 \mu g \ of \ Lyophilised \ peptide, \ formulation \ De-ionized \ (DI) \ water.$

Unit Size 100 μα

Storage

Shipped at ambient temperature, store at -20°C. Instructions

Accession

NP_780512.2; NP_958797.2; NP_958798.2 ID

Peptide with sequence C-DQSKKGKGKPQN, from the internal region of the protein sequence according to Immunogen

NP 780512.2; NP 958797.2; NP 958798.2.

The peptide was used in the production of EB07782 - Goat Anti-Sall4 (mouse) Antibody. It is currently untested **Product**

Comments in blocking applications.

Peptide

C-DQSKKGKGKPQN Sequence

Shipping

Instructions Ambient

Product Grade

https://prod-vector-labs-pimcore-assets.s3.us-east-1.amazonaws.com/assets/products/image/elite_medium.png