

Email: customerservice@vectorlabs.com

Telephone: (650) 697-3600

EBP07813 - IMMUNIZING PEPTIDE

SKU: EBP07813

SPECIFICATIONS

 $\textbf{Formulation} \ 100 \mu g \ of \ Lyophilised \ peptide, \ formulation \ De-ionized \ (DI) \ water.$

Unit Size 100 μα

Storage

Shipped at ambient temperature, store at -20°C. Instructions

Accession

NP_001036214.1; NP_001036215.1; NP_001036216.1; NP_006594.3 ID

Peptide with sequence C-REPKRLRPEDS, from the internal region of the protein sequence according to Immunogen

NP 001036214.1; NP 001036215.1; NP 001036216.1; NP 006594.3.

The peptide was used in the production of EB07813 - Goat Anti-STAG2 Antibody. It is currently untested in **Product**

Comments blocking applications.

Peptide

C-REPKRLRPEDS Sequence

Shipping

Instructions Ambient

Product Grade

https://prod-vector-labs-pimcore-assets.s3.us-east-1.amazonaws.com/assets/products/image/elite_medium.png