

# EBP07813 - IMMUNIZING PEPTIDE

SKU: EBP07813



---

## SPECIFICATIONS

<b>Formulation</b>	100µg of Lyophilised peptide, formulation De-ionized (DI) water.
<b>Unit Size</b>	100 µg
<b>Storage Instructions</b>	Shipped at ambient temperature, store at -20°C.
<b>Accession ID</b>	NP_001036214.1; NP_001036215.1 ; NP_001036216.1 ; NP_006594.3
<b>Immunogen</b>	Peptide with sequence C-REPKRLRPEDS, from the internal region of the protein sequence according to NP_001036214.1; NP_001036215.1 ; NP_001036216.1 ; NP_006594.3.
<b>Product Comments</b>	The peptide was used in the production of EB07813 - Goat Anti-STAG2 Antibody. It is currently untested in blocking applications.
<b>Peptide Sequence</b>	C-REPKRLRPEDS
<b>Shipping Instructions</b>	Ambient
<b>Product Grade</b>	<a href="https://prod-vector-labs-pimcore-assets.s3.us-east-1.amazonaws.com/assets/products/image/elite_medium.png">https://prod-vector-labs-pimcore-assets.s3.us-east-1.amazonaws.com/assets/products/image/elite_medium.png</a>

## DOCUMENTS

- [Data Sheet](#)