

Email: customerservice@vectorlabs.com

Telephone: (650) 697-3600

EBP07819 - IMMUNIZING PEPTIDE

SKU: EBP07819

SPECIFICATIONS

Formulation 100µg of Lyophilised peptide, formulation De-ionized (DI) water.

Unit Size 100 μg

Storage

Instructions Shipped at ambient temperature, store at -20°C.

Accession

NP_766100.1; NP_084161.1 ID

Immunogen Peptide with sequence C-KSGRTTEPHTN, from the C Terminus of the protein sequence according to NP_766100.1; NP_084161.1.

The peptide was used in the production of EB07819 - Goat Anti-Gabpb2 (mouse) Antibody. It is currently untested Product

Comments in blocking applications.

Peptide

C-KSGRTTEPHTN Sequence

Shipping

Ambient Instructions

Product

https://prod-vector-labs-pimcore-assets.s3.us-east-1.amazonaws.com/assets/products/image/aspiring_medium.png Grade