

EBP07846 - IMMUNIZING PEPTIDE

SKU: EBP07846



SPECIFICATIONS

Formulation	100µg of Lyophilised peptide, formulation De-ionized (DI) water.
Unit Size	100 µg
Storage Instructions	Shipped at ambient temperature, store at -20°C.
Accession ID	NP_058590.1; NP_113833.1
Immunogen	Peptide with sequence CDWGKTIQSYGE, from the internal region of the protein sequence according to NP_058590.1; NP_113833.1.
Product Comments	The peptide was used in the production of EBP07846 - Goat Anti-Ramp1 (mouse) (aa57-68) Antibody. It is currently untested in blocking applications.
Peptide Sequence	CDWGKTIQSYGE
Shipping Instructions	Ambient
Product Grade	https://prod-vector-labs-pimcore-assets.s3.us-east-1.amazonaws.com/assets/products/image/aspiring_medium.png