

Email: customerservice@vectorlabs.com

Telephone: (650) 697-3600

EBP07846 - IMMUNIZING PEPTIDE

SKU: EBP07846

SPECIFICATIONS

Formulation 100µg of Lyophilised peptide, formulation De-ionized (DI) water.

Unit Size 100 μg

Storage

Instructions Shipped at ambient temperature, store at -20°C.

Accession

NP_058590.1; NP_113833.1 ID

 $\textbf{Immunogen} \begin{tabular}{ll} \textbf{Peptide with sequence CDWGKTIQSYGE, from the internal region of the protein sequence according to NP_058590.1; NP_113833.1. \end{tabular}$

The peptide was used in the production of EB07846 - Goat Anti-Ramp1 (mouse) (aa57-68) Antibody. It is currently Product

Comments untested in blocking applications.

Peptide

CDWGKTIQSYGE Sequence

Shipping

Ambient Instructions

Product

https://prod-vector-labs-pimcore-assets.s3.us-east-1.amazonaws.com/assets/products/image/aspiring_medium.png Grade