

EBP07847 - IMMUNIZING PEPTIDE

SKU: EBP07847



SPECIFICATIONS

Formulation	100µg of Lyophilised peptide, formulation De-ionized (DI) water.
Unit Size	100 µg
Storage	
Instructions	Shipped at ambient temperature, store at -20°C.
Accession ID	NP_006051.1; NP_001207694.1; NP_001207695.1; NP_001207696.1; NP_001207697.1; NP_001207698.1; NP_001207699.1; NP_001207700.1 and more
Immunogen	Peptide with sequence C-DLLRAASENSQDALR, from the internal region of the protein sequence according to NP_006051.1; NP_001207694.1; NP_001207695.1; NP_001207696.1; NP_001207697.1; NP_001207698.1; NP_001207699.1; NP_001207700.1 and more.
Product Comments	The peptide was used in the production of EBP07847 - Goat Anti-IKZF1 / IKAROS Antibody. It is currently untested in blocking applications.
Peptide Sequence	C-DLLRAASENSQDALR
Shipping Instructions	Ambient
Product Grade	https://prod-vector-labs-pimcore-assets.s3.us-east-1.amazonaws.com/assets/products/image/elite_medium.png

DOCUMENTS

- [Data Sheet](#)