

Email: customerservice@vectorlabs.com

Telephone: (650) 697-3600

EBP07847 - IMMUNIZING PEPTIDE

SKU: EBP07847

SPECIFICATIONS

 $\textbf{Formulation} \ 100 \mu g \ of \ Lyophilised \ peptide, \ formulation \ De-ionized \ (DI) \ water.$

Unit Size 100 μα

Shipped at ambient temperature, store at -20°C.

NP_006051.1; NP_001207694.1; NP_001207695.1; NP_001207696.1; NP_001207697.1; NP_001207698.1; Accession

NP_001207699.1; NP_001207700.1 and more ID

Peptide with sequence C-DLLRAASENSQDALR, from the internal region of the protein sequence according to

Immunogen NP 006051.1; NP 001207694.1; NP 001207695.1; NP 001207696.1; NP 001207697.1; NP 001207698.1;

NP_001207699.1; NP_001207700.1 and more.

The peptide was used in the production of EB07847 - Goat Anti-IKZF1 / IKAROS Antibody. It is currently **Product**

Comments untested in blocking applications.

Peptide

C-DLLRAASENSQDALR Sequence

Shipping

Instructions Ambient

Product

https://prod-vector-labs-pimcore-assets.s3.us-east-1.amazonaws.com/assets/products/image/elite_medium.png Grade