

Email: customerservice@vectorlabs.com

Telephone: (650) 697-3600

EBP07850 - IMMUNIZING PEPTIDE

SKU: EBP07850

SPECIFICATIONS

Formulation 100µg of Lyophilised peptide, formulation De-ionized (DI) water.

Unit Size 100 μg

Storage

Instructions Shipped at ambient temperature, store at -20°C.

Accession

NP_002362.1; NP_071884.1 ID

The peptide was used in the production of EB07850 - Goat Anti-T-cell differentiation protein Antibody. It is Product

Comments currently untested in blocking applications.

Peptide

C-QDGFTYRHYHEN Sequence

Shipping

Ambient Instructions

Product

https://prod-vector-labs-pimcore-assets.s3.us-east-1.amazonaws.com/assets/products/image/aspiring_medium.png Grade