

EBP07856 - IMMUNIZING PEPTIDE

SKU: EBP07856



SPECIFICATIONS

| | |
|------------------------------|---|
| Formulation | 100µg of Lyophilised peptide, formulation De-ionized (DI) water. |
| Unit Size | 100 µg |
| Storage Instructions | Shipped at ambient temperature, store at -20°C. |
| Accession ID | NP_037393.1; NP_001317680.1; NP_001317681.1; NP_001317682.1; NP_001341756.1; |
| Immunogen | Peptide with sequence C-DGLFDDSEDESDK, from the internal region of the protein sequence according to NP_037393.1; NP_001317680.1; NP_001317681.1; NP_001317682.1; NP_001341756.1;. |
| Product Comments | The peptide was used in the production of EB07856 - Goat Anti-PGC1A/PPARGC1A Antibody. It is currently untested in blocking applications. |
| Peptide Sequence | C-DGLFDDSEDESDK |
| Shipping Instructions | Ambient |
| Product Grade | https://prod-vector-labs-pimcore-assets.s3.us-east-1.amazonaws.com/assets/products/image/elite_plus_medium.png |

DOCUMENTS

- [Data Sheet](#)