

# EBP07865 - IMMUNIZING PEPTIDE

**SKU:** EBP07865



---

## SPECIFICATIONS

|                              |                                                                                                                                                                                                                                   |
|------------------------------|-----------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| <b>Formulation</b>           | 100µg of Lyophilised peptide, formulation De-ionized (DI) water.                                                                                                                                                                  |
| <b>Unit Size</b>             | 100 µg                                                                                                                                                                                                                            |
| <b>Storage</b>               |                                                                                                                                                                                                                                   |
| <b>Instructions</b>          | Shipped at ambient temperature, store at -20°C.                                                                                                                                                                                   |
| <b>Accession ID</b>          | NP_444253.3; NP_444254.3 ; NP_444255.3; NP_444256.3                                                                                                                                                                               |
| <b>Immunogen</b>             | Peptide with sequence C-KDTKNMEAKKLSKD, from the internal region of the protein sequence according to NP_444253.3; NP_444254.3 ; NP_444255.3; NP_444256.3.                                                                        |
| <b>Product</b>               | The peptide was used in the production of EBP07865 - Goat Anti-MYLK Antibody. It is currently untested in                                                                                                                         |
| <b>Comments</b>              | blocking applications.                                                                                                                                                                                                            |
| <b>Peptide Sequence</b>      | C-KDTKNMEAKKLSKD                                                                                                                                                                                                                  |
| <b>Shipping Instructions</b> | Ambient                                                                                                                                                                                                                           |
| <b>Product Grade</b>         | <a href="https://prod-vector-labs-pimcore-assets.s3.us-east-1.amazonaws.com/assets/products/image/elite_medium.png">https://prod-vector-labs-pimcore-assets.s3.us-east-1.amazonaws.com/assets/products/image/elite_medium.png</a> |

## DOCUMENTS

- [Data Sheet](#)