

Email: customerservice@vectorlabs.com

Telephone: (650) 697-3600

EBP07870 - IMMUNIZING PEPTIDE

SKU: EBP07870

SPECIFICATIONS

 $\textbf{Formulation} \ 100 \mu g \ of \ Lyophilised \ peptide, \ formulation \ De-ionized \ (DI) \ water.$

Unit Size 100 μα

Shipped at ambient temperature, store at -20°C.

NP_001123645.1; NP_005366.2; NP_001123644.1; NP_001155128.1; NP_001155129.1; NP_001155130.1; Accession

ID NP_001155131.1

Peptide with sequence C-QRHYNDEDPEKEKR, from the internal region of the protein sequence according to

Immunogen NP 001123645.1; NP 005366.2; NP 001123644.1; NP 001155128.1; NP 001155129.1; NP 001155130.1;

NP 001155131.1.

The peptide was used in the production of EB07870 - Goat Anti-c-Myb (aa281-294) Antibody. It is currently **Product**

Comments untested in blocking applications.

Peptide

C-QRHYNDEDPEKEKR Sequence

Shipping

Instructions Ambient

Product

https://prod-vector-labs-pimcore-assets.s3.us-east-1.amazonaws.com/assets/products/image/elite_medium.png Grade