



Email: <a href="mailto:customerservice@vectorlabs.com">customerservice@vectorlabs.com</a>

Telephone: (650) 697-3600

## **EBP07872 - IMMUNIZING PEPTIDE**

**SKU:** EBP07872

## **SPECIFICATIONS**

Formulation 100µg of Lyophilised peptide, formulation De-ionized (DI) water.

**Unit Size** 100 μg

Storage

Instructions Shipped at ambient temperature, store at -20°C.

**Accession** 

NP\_032343.2 ID

Immunogen Peptide with sequence SDKTSTKESCQ, from the internal region of the protein sequence according to NP\_032343.2.

Product The peptide was used in the production of EB07872 - Goat Anti-hyaluronidase 1 (mouse, rat) Antibody. It is

**Comments** currently untested in blocking applications.

**Peptide** 

**SDKTSTKESCQ** Sequence

**Shipping** 

Ambient Instructions

**Product** 

https://prod-vector-labs-pimcore-assets.s3.us-east-1.amazonaws.com/assets/products/image/aspiring medium.png Grade