

Email: <a href="mailto:customerservice@vectorlabs.com">customerservice@vectorlabs.com</a>

Telephone: (650) 697-3600

## **EBP07883 - IMMUNIZING PEPTIDE**

**SKU:** EBP07883

## **SPECIFICATIONS**

 $\textbf{Formulation} \ 100 \mu g \ of \ Lyophilised \ peptide, \ formulation \ De-ionized \ (DI) \ water.$ 

**Unit Size** 100 μα

Shipped at ambient temperature, store at -20°C.

Accession

 $NP\_002533.1; \ NP\_058212.1; \ NP\_058213.1; \ NP\_058214.1; \ NP\_058434.1; \ NP\_058436.1; \ NP\_058437.1; \\$ ID NP\_058438.1

Peptide with sequence C-DLRQSRHAQEPPAK, from the C Terminus of the protein sequence according to

Immunogen NP 002533.1; NP 058212.1; NP 058213.1; NP 058214.1; NP 058434.1; NP 058436.1; NP 058437.1;

The peptide was used in the production of EB07883 - Goat Anti-OGG1 Antibody. It is currently untested in **Product** 

Comments blocking applications.

**Peptide** 

C-DLRQSRHAQEPPAK Sequence

Shipping

Instructions Ambient

**Product** 

https://prod-vector-labs-pimcore-assets.s3.us-east-1.amazonaws.com/assets/products/image/elite\_medium.png Grade