

EBP07883 - IMMUNIZING PEPTIDE

SKU: EBP07883



SPECIFICATIONS

Formulation 100µg of Lyophilised peptide, formulation De-ionized (DI) water.

Unit Size 100 µg

Storage Instructions Shipped at ambient temperature, store at -20°C.

Accession ID NP_002533.1; NP_058212.1 ; NP_058213.1; NP_058214.1; NP_058434.1; NP_058436.1; NP_058437.1; NP_058438.1

Immunogen Peptide with sequence C-DLRQSRHAQEPPAK, from the C Terminus of the protein sequence according to NP_002533.1; NP_058212.1 ; NP_058213.1; NP_058214.1; NP_058434.1; NP_058436.1; NP_058437.1; NP_058438.1.

Product Comments The peptide was used in the production of EBP07883 - Goat Anti-OGG1 Antibody. It is currently untested in blocking applications.

Peptide Sequence C-DLRQSRHAQEPPAK

Shipping Instructions Ambient

Product Grade https://prod-vector-labs-pimcore-assets.s3.us-east-1.amazonaws.com/assets/products/image/elite_medium.png