

EBP07896 - IMMUNIZING PEPTIDE

SKU: EBP07896



SPECIFICATIONS

Formulation	100µg of Lyophilised peptide, formulation De-ionized (DI) water.
Unit Size	100 µg
Storage Instructions	Shipped at ambient temperature, store at -20°C.
Accession ID	NP_003378.3; NP_001070737.1; NP_001362764.1; NP_001362768.1
Immunogen	Peptide with sequence C-KLARNESRSGSNRRE, from the C Terminus of the protein sequence according to NP_003378.3; NP_001070737.1; NP_001362764.1; NP_001362768.1.
Product Comments	The peptide was used in the production of EBP07896 - Goat Anti-WIPF1 Antibody. It is currently untested in blocking applications.
Peptide Sequence	C-KLARNESRSGSNRRE
Shipping Instructions	Ambient
Product Grade	https://prod-vector-labs-pimcore-assets.s3.us-east-1.amazonaws.com/assets/products/image/elite_medium.png

DOCUMENTS

- [Data Sheet](#)