

EBP07898-B - IMMUNIZING PEPTIDE

SKU: EBP07898-B



SPECIFICATIONS

Formulation	100µg of Lyophilised peptide, formulation De-ionized (DI) water.
Unit Size	100 µg
Storage Instructions	Shipped at ambient temperature, store at -20°C.
Accession ID	NP_004960.2; NP_001159418.1; NP_001309722.1; NP_001309723.1; NP_001309724.1; NP_001309726.1
Immunogen	Peptide with sequence CQQYNFDRDNTE., from the internal region of the protein sequence according to NP_004960.2; NP_001159418.1; NP_001309722.1; NP_001309723.1; NP_001309724.1; NP_001309726.1.
Product Comments	The peptide was used in the production of EBP07898-B - Goat Anti-Insulysin / Insulinase, Biotinylated Antibody. It is currently untested in blocking applications.
Peptide Sequence	CQQYNFDRDNTE.
Shipping Instructions	Ambient
Product Grade	https://prod-vector-labs-pimcore-assets.s3.us-east-1.amazonaws.com/assets/products/image/elite_medium.png

DOCUMENTS

- [Data Sheet](#)