

Email: <a href="mailto:customerservice@vectorlabs.com">customerservice@vectorlabs.com</a>

Telephone: (650) 697-3600

## **EBP07903 - IMMUNIZING PEPTIDE**

**SKU:** EBP07903

## **SPECIFICATIONS**

Formulation 100µg of Lyophilised peptide, formulation De-ionized (DI) water.

**Unit Size** 100 μg

Storage

Instructions Shipped at ambient temperature, store at -20°C.

Accession

NP\_996800.1; NP\_996798.1; NP\_055195.1 ID

The peptide was used in the production of EB07903 - Goat Anti-NGFRAP1 Antibody. It is currently untested in Product

**Comments** blocking applications. **Peptide** 

Sequence

**CLRILMGELSNHHDH** 

Shipping

Ambient Instructions

**Product** 

https://prod-vector-labs-pimcore-assets.s3.us-east-1.amazonaws.com/assets/products/image/aspiring\_medium.png Grade