

EBP07950 - IMMUNIZING PEPTIDE

SKU: EBP07950



SPECIFICATIONS

Formulation	100µg of Lyophilised peptide, formulation De-ionized (DI) water.
Unit Size	100 µg
Storage Instructions	Shipped at ambient temperature, store at -20°C.
Accession ID	NP_789829.1
Immunogen	Peptide with sequence C-NPNSHSETGSKAGDK, from the C Terminus of the protein sequence according to NP_789829.1.
Product Comments	The peptide was used in the production of EBP07950 - Goat Anti-Tfap2a (zebrafish) Antibody. It is currently untested in blocking applications.
Peptide Sequence	C-NPNSHSETGSKAGDK
Shipping Instructions	Ambient
Product Grade	https://prod-vector-labs-pimcore-assets.s3.us-east-1.amazonaws.com/assets/products/image/aspiring_medium.png

DOCUMENTS

- [Data Sheet](#)