

Email: <a href="mailto:customerservice@vectorlabs.com">customerservice@vectorlabs.com</a>

Telephone: (650) 697-3600

## **EBP07957 - IMMUNIZING PEPTIDE**

**SKU:** EBP07957

## **SPECIFICATIONS**

 $\textbf{Formulation} \ 100 \mu g \ of \ Lyophilised \ peptide, \ formulation \ De-ionized \ (DI) \ water.$ 

**Unit Size** 100 μg

Storage

Instructions Shipped at ambient temperature, store at -20°C.

**Accession** 

NP\_006152.2 ID

Immunogen Peptide with sequence KDLLHTTPCFIPYH, from the C Terminus of the protein sequence according to NP\_006152.2.

Product The peptide was used in the production of EB07957 - Goat Anti-Neurogenin 1 Antibody. It is currently untested in

**Comments** blocking applications.

**Peptide** 

Sequence

**KDLLHTTPCFIPYH** 

**Shipping** Ambient Instructions

**Product** 

Grade

https://prod-vector-labs-pimcore-assets.s3.us-east-1.amazonaws.com/assets/products/image/aspiring\_medium.png