

Email: customerservice@vectorlabs.com

Telephone: (650) 697-3600

EBP07976 - IMMUNIZING PEPTIDE

SKU: EBP07976

SPECIFICATIONS

 $\textbf{Formulation} \ 100 \mu g \ of \ Lyophilised \ peptide, \ formulation \ De-ionized \ (DI) \ water.$

Unit Size 100 μα

Shipped at ambient temperature, store at -20°C.

NP_647593.1; NP_647594.1; NP_647595.1; NP_647596.1; NP_647597.1; NP_647598.1; NP_647599.1; Accession

NP_004296.1; NP_647600.1; NP_647601.1 ID

Peptide with sequence C-KESDWNQHKELEK, from the internal region (near the C Terminus) of the protein

Immunogen sequence according to NP 647593.1; NP 647594.1; NP 647595.1; NP 647596.1; NP 647597.1; NP 647598.1;

NP_647599.1; NP_004296.1; NP_647600.1; NP_647601.1.

The peptide was used in the production of EB07976 - Goat Anti-Amphiphysin2 / BIN1 Antibody. It is currently **Product**

Comments untested in blocking applications.

Peptide

C-KESDWNQHKELEK Sequence

Shipping

Instructions Ambient

Product

https://prod-vector-labs-pimcore-assets.s3.us-east-1.amazonaws.com/assets/products/image/elite_medium.png Grade