

EBP07976 - IMMUNIZING PEPTIDE

SKU: EBP07976

SPECIFICATIONS

Formulation	100µg of Lyophilised peptide, formulation De-ionized (DI) water.
Unit Size	100 µg
Storage	
Instructions	Shipped at ambient temperature, store at -20°C.
Accession ID	NP_647593.1; NP_647594.1; NP_647595.1; NP_647596.1; NP_647597.1; NP_647598.1; NP_647599.1; NP_004296.1; NP_647600.1; NP_647601.1
Immunogen	Peptide with sequence C-KESDWNQHKELEK, from the internal region (near the C Terminus) of the protein sequence according to NP_647593.1; NP_647594.1; NP_647595.1; NP_647596.1; NP_647597.1; NP_647598.1; NP_647599.1; NP_004296.1; NP_647600.1; NP_647601.1.
Product Comments	The peptide was used in the production of EBP07976 - Goat Anti-Amphiphysin2 / BIN1 Antibody. It is currently untested in blocking applications.
Peptide Sequence	C-KESDWNQHKELEK
Shipping Instructions	Ambient
Product Grade	https://prod-vector-labs-pimcore-assets.s3.us-east-1.amazonaws.com/assets/products/image/elite_medium.png