

# EBP07979 - IMMUNIZING PEPTIDE

**SKU:** EBP07979



## SPECIFICATIONS

<b>Formulation</b>	100µg of Lyophilised peptide, formulation De-ionized (DI) water.
<b>Unit Size</b>	100 µg
<b>Storage</b>	Shipped at ambient temperature, store at -20°C.
<b>Instructions</b>	
<b>Accession ID</b>	NP_001196.1; NP_053581.1; NP_053582.1; NP_053583.1
<b>Immunogen</b>	Peptide with sequence C-DLLRQRKTTKESLAQ, from the internal region of the protein sequence according to NP_001196.1; NP_053581.1; NP_053582.1; NP_053583.1.
<b>Product</b>	The peptide was used in the production of EBP07979 - Goat Anti-BNIP1 Antibody. It is currently untested in blocking applications.
<b>Comments</b>	
<b>Peptide Sequence</b>	C-DLLRQRKTTKESLAQ
<b>Shipping Instructions</b>	Ambient
<b>Product Grade</b>	<a href="https://prod-vector-labs-pimcore-assets.s3.us-east-1.amazonaws.com/assets/products/image/elite_medium.png">https://prod-vector-labs-pimcore-assets.s3.us-east-1.amazonaws.com/assets/products/image/elite_medium.png</a>

## DOCUMENTS

- [Data Sheet](#)