

Email: <a href="mailto:customerservice@vectorlabs.com">customerservice@vectorlabs.com</a>

Telephone: (650) 697-3600

## **EBP08000 - IMMUNIZING PEPTIDE**

**SKU:** EBP08000

## **SPECIFICATIONS**

 $\textbf{Formulation} \ 100 \mu g \ of \ Lyophilised \ peptide, \ formulation \ De-ionized \ (DI) \ water.$ 

**Unit Size** 100 μα

Storage

Shipped at ambient temperature, store at -20°C. Instructions

**Accession** 

NP\_000821.1; NP\_783300.1 ID

Peptide with sequence QCKQTHPTNSTS, from the internal region (near the C Terminus) of the protein Immunogen

sequence according to NP 000821.1; NP 783300.1.

The peptide was used in the production of EB08000 - Goat Anti-GLuR5 / GRIK1 Antibody. It is currently **Product** 

**Comments** untested in blocking applications.

**Peptide** 

**QCKQTHPTNSTS** Sequence

Shipping

Instructions Ambient

**Product** Grade

https://prod-vector-labs-pimcore-assets.s3.us-east-1.amazonaws.com/assets/products/image/elite\_medium.png